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GENOMIC FINDINGS	THERAPIES WITH CLINICAL BENEFIT (IN PATIENT'S TUMOR TYPE)
<i>BRAF</i> - V600E	Dabrafenib + Trametinib <span style="border: 1px solid black; padding: 0 2px;">2A</span>

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GENE ALTERATIONS	VAF %	THERAPIES APPROVED IN THE EU (IN PATIENT'S TUMOR TYPE)	THERAPIES APPROVED IN THE EU (IN OTHER TUMOR TYPE)
<i>BRCA1</i> - complex rearrangement	19.4%	Olaparib <span style="border: 1px solid black; padding: 0 2px;">1</span>	Niraparib
<i>ET1660fs*17</i>	20.2%	Talazoparib <span style="border: 1px solid black; padding: 0 2px;">1</span>	Rucaparib

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**VARIANTS TO CONSIDER FOR FOLLOW-UP GERMLINE TESTING IN SELECT CANCER SUSCEPTIBILITY GENES**

*Findings below have been previously reported as pathogenic germline in the ClinVar genomic data and were detected at an allele frequency of >30%. See appendix for details.*

**BRCA2 - S1982fs\*22** ..... p. 5

This report does not indicate whether variants listed above are germline or somatic in this patient. In the appropriate clinical context, follow-up germline testing would be needed to determine whether a finding is germline or somatic.

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ESMO: European Society for Medical Oncology; PMWG: Precision Medicine Working Group.

## Resistenzen

GENOMIC FINDINGS	THERAPIES WITH CLINICAL BENEFIT (IN PATIENT'S TUMOR TYPE)	THERAPIES WITH CLINICAL BENEFIT (IN OTHER TUMOR TYPE)
<b>ERBB2 - L755S, amplification</b>	Ado-trastuzumab emtansine <span style="border: 1px solid gray; padding: 2px;">1</span>	Afatinib
	Pertuzumab <span style="border: 1px solid gray; padding: 2px;">1</span>	Dacomitinib
	Neratinib <span style="border: 1px solid gray; padding: 2px;">2A</span>	
	Fam-trastuzumab deruxtecan <span style="border: 1px solid gray; padding: 2px;">2A</span>	
	Lapatinib <span style="border: 1px solid gray; padding: 2px;">?</span>	
	Trastuzumab <span style="border: 1px solid gray; padding: 2px;">x</span>	
	10 Trials see p.26	

? Limited evidence showing variant(s) in this sample may confer resistance to this therapy

x Extensive evidence showing variant(s) in this sample may confer resistance to this therapy

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GENOMIC SIGNATURES	THERAPIES APPROVED IN THE EU (IN PATIENT'S TUMOR TYPE)	THERAPIES APPROVED IN THE EU (IN OTHER TUMOR TYPE)
Loss of Heterozygosity score - 27.0 %	Niraparib [2A]	Talazoparib
	Olaparib [2A]	
	Rucaparib [2A]	

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Tumor Fraction - 22%	Tumor fraction is an estimate of the percentage of circulating-tumor DNA (ctDNA) present in a cell-free DNA (cfDNA) sample based on observed aneuploid instability.
----------------------	---

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
Blood Tumor Mutational Burden - 6 Muts/Mb	No therapies or clinical trials. See Genomic Signatures section
---	---

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**Microsatellite status - MSI-High Not Detected**

**MSI-High not detected. No evidence of microsatellite instability in this sample (see Appendix section).**

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### Referenzen:

1. <https://www.foundationmedicine.at/> (Accessed on April 2021)
2. Data on File: FoundationOne® CDx Sample Report 2021 updates, 2, Mandelker D et al. Ann Oncol. 2019; 30: 1221-1231.
3. Bronkhorst AJ, et al. Biomol Detect Quantif (2019) PMID: 30923679
4. Raja R, et al. Clin. Cancer Res. (2018) PMID: 30093454
5. Hrebien S, et al. Ann. Oncol. (2019) PMID: 30860573
6. Choudhury AD, et al. JCI Insight (2018) PMID: 30385733
7. Goodall J, et al. Cancer Discov (2017) PMID: 28450425
8. Goldberg SB, et al. Clin. Cancer Res. (2018) PMID: 29330207
9. Bettegowda C, et al. Sci Transl Med (2014) PMID: 24553385
10. Hall MJ, et al. Evaluation of microsatellite instability (MSI) status 11,573 diverse solid tumors using comprehensive genomic profiling (CGP). Poster session presented at: 2016 ASCO Annual Meeting; 2016 June 3-7; Chicago, IL.
11. Trabucco SE, Gowen K, Maund SL, et al. A Novel Next-Generation Sequencing Approach to Detecting Microsatellite Instability and Pan-Tumor Characterization of 1000 Microsatellite Instability-High Cases in 67,000 Patient Samples. J Mol Diagn. 2019;21(6):1053-1066. doi:10.1016/j.jmoldx.2019.06.011
12. Middha S, Zhang L, Nafa K, et al. Reliable Pan-Cancer Microsatellite Instability Assessment by Using Targeted Next-Generation Sequencing Data. JCO Precis Oncol. 2017;2017.10.1200/PO.17.00084. doi:10.1200/PO.17.00084
13. Trabucco SE, Gowen K, Maund SL, et al. A Novel Next-Generation Sequencing Approach to Detecting Microsatellite Instability and Pan-Tumor Characterization of 1000 Microsatellite Instability-High Cases in 67,000 Patient Samples. J Mol Diagn. 2019;21(6):1053-1066. doi:10.1016/j.jmoldx.2019.06.011

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© 2021 Foundation Medicine, Inc. Foundation Medicine®, FoundationOne® CDx, FoundationOne® LIQUID CDx und FoundationOne® Heme sind eingetragene Warenzeichen. Roche ist der lizenzierte Anbieter von Foundation Medicine Produkten außerhalb der Vereinigten Staaten von Amerika.